


# Bench Alignment


Lock the knee from the app.

- 1 Reference Line:**  
Create a vertical Reference Line (A) using a plumb.
- 2 Foot Placement:**  
Keep the shoe on. Keep Mid-foot approximately  $30 \pm 10$  mm anterior (I) to the Reference Line (A). Place the foot with a  $5-7^\circ$  toe-out angle.
- 3 Knee Placement:**  
Place the knee with approximately  $5^\circ$  external rotation. Adjust the fore-aft set screws of the foot to move the Knee Axis (H) within  $\pm 5$  mm w.r.t. the Reference Line (A).
- 4 Mark the Socket:**  
Mark the Proximal Center (B), Distal Center (E), and the Socket Centerline (C).
- 5 Socket Alignment:**  
Let the Reference Line (A) pass the Proximal Center (B) and set the Socket Flexion Angle (D) to  $3-5^\circ$ . Take individual conditions into account, such as hip flexion contractures.
- 6 Socket Connection:**  
Align and fixate the Socket Adapter (F) to the socket and connect the knee to the socket via the Socket Adapter (F).



# Tips for alignment

Neutral & dynamic alignment is recommended.



Neutral & dynamic    Overly-stable

**Overly-stable alignment could decrease functionality.**

It could potentially make it difficult to trigger stair ascent and descent, stand-to-sit, ramp descent and walk or trigger swing.

A. Reference line, B. Proximal Center, C. Socket Centerline, D. Socket Flexion Angle, E. Distal Center, F. Socket Adapter, G. Tube Adapter, H. Knee Axis, I. Distance between the reference line and mid-foot.

## Static Alignment

Lock the knee from the app.


Ask the patient to put on the prosthesis.

### Without the 3D L.A.S.A.R. Posture tool

- 1 **Stand Straight**  
Ask the patient to stand straight and distribute weight evenly on both feet.
- 2 **Foot Adjustment**  
Adjust solely the plantar- and dorsiflexion of the foot so that the pressure on the front or rear side of the socket are even.

### With the 3D L.A.S.A.R. Posture tool

- 1 **Load Line**  
Instruct the wearer to step onto the force plate of the 3D L.A.S.A.R. Posture tool with the prosthetic foot and place the other foot on a height compensation panel. Instruct the wearer to distribute weight evenly on both feet.  
The Load Line will be visible on the side.
- 2 **Foot Adjustment**  
Adjust only the plantar or dorsiflexion of the foot so that the Knee Axis mark is  $30 \pm 10$  mm posterior to the Load Line.



J. Loading Line, K. Distance between the Knee Axis and Loading Line, L. Adapter for the Foot.

## Dynamic Alignment

Set the knee to Start mode.

- 1 **Sagittal Plane**  
Check stride symmetry and stance stability, adjust knee flexion or foot flexion when necessary.
- 2 **Frontal Plane**  
Observe the lateral motion and stability of trunk. Check pressure on medial-lateral side of the socket with the wearer. Adjust ab/adduction and position of the socket adapter when necessary. Adjust ankle inverse of eversion when necessary.
- 3 **Transversal Plane**  
Adjust knee rotation when excessive leg rotation is observed during swing. Adjust foot rotation when leg rotation is observed e.s.p. during stance.