

How to Adjust the Settings in the Extra Tab?



- 1 Reset Defaults**
Reset the parameters in this Extra Tab to their factory default.
- 2 Activity Level**
Adjust the activity level according to the user profile. The higher the value, the easier it is to achieve a higher walking speed. 3 or 4 are suitable for most users.
- 3 Extension Speed**
Adjusting the extension speed impacts how early swing extension is triggered. The higher the value, the sooner extending starts during walking.
- 4 Force Sensor Sensitivity**
The higher the value, the easier it is to trigger swing.
- 5 Flexion Amplitude**
Flexion amplitude is the maximum flexion angle that Intuy Knee will reach while walking. The bigger the value, the lighter it feels to bring the prosthetic leg forward. 65-70 degrees is recommended.
- 6 Pre-flex Angle**
Pre-flex angle is how far the knee is extending during terminal swing. The higher the value, the more stance yielding.
- 7 Stance Resistance**
Stance resistance affects sitting down and descending the stairs/ramps. The higher the value, the more resistance.
- 8 Sit-to-Stand Power Level**
The higher the value, the more support user gets while standing up.
- 9 Stair Ascent Power Level**
The higher the value, the more support user gets while climbing the stairs.
- 10 Enable Stair Mode**
You are able to enable and disable the stair mode if necessary, by checking the box.
- 11 Save Settings**
Click the checkmark after making any changes to save them.

Note: When adjusting these parameters, make sure the user is aware of the change. Perform the adjustment in a safe environment, such as on level ground not on stairs, while standing still and not while walking.